


Guidelines for FibroTrack installation software for windows users

To install FibroTrack on windows, open


https://drive.google.com/file/d/15VkjJnp3pEWw3y8nUWc_g3SFsi_UiOKi/view?usp=sharing

Click on the download button (Step1) and then Download anyway (Step 2). Finally, unzip the file by pressing “extract all” or “extract here” button (Step3).


Step 1



Step 2



Step 3

















Step 4

Installation steps:


- Open the extracted “FibroTrack installation” folder and then open the “FibroTrack - Software Icons” folder.
- Find "Install FibroTrack Graphical User Interface.exe" (A)
- Double-click to run. You will have a pop up message: “Windows Protected your PC” click on “More info” (B).
- Click "Run anyway" (C).

Note: Do not delete any files from this folder - all contents are required for proper installation!


A

Name	Date	Type
 best.pt	24/10/2024 22:01	PT File
 FibroTrack ICON	30/10/2024 7:46	ICO File
 FibroTrack_Final_Version_Source_Code	06/01/2025 15:36	Application
 FibroTrack_Final_Version_Source_Code.py	01/01/2025 16:10	PY File
 FibroTrackIcon	27/08/2024 21:36	PNG File
 Group 1	27/08/2024 19:11	PNG File
 Group 2	27/08/2024 19:10	PNG File
 Install FibroTrack Graphical User Interface	06/01/2025 15:41	Application
 Install FibroTrack Graphical User Interface-1a.bin	06/01/2025 15:41	BIN File
 Odeh et al. Analyze Fibrosis	26/08/2024 14:53	PNG File
 Odeh et al. Define Pixel	26/08/2024 12:27	PNG File
 Odeh et al. Muscle Segmentation	26/08/2024 14:33	PNG File
 Splash Screen	27/08/2024 18:43	TIF File
 Splash_Screen_FibroTrack	01/09/2024 20:03	TIF File

B






C




Step 5

Install the software by accepting the license.





The splash screen appears with a black console. It may take a few minutes to launch the software while importing all the modules for deep learning.



Finally, a black console will appear with the following message.


```
Initializing the application...
please wait...

Welcome To

Imports completed.
Please wait, you're almost there...
Launching the FibroTrack application...
```

Note: Don't close this console while you're using the software, otherwise you will close the software.

In the meanwhile, the application will be popped up and cover the whole screen of your computer as follows:



Now, you are ready to use FibroTrack.