




Guidelines for FibroTrack installation software for windows users



To install FibroTrack on windows, open <https://github.com/Anas-Odeh/FibroTrack>. Click on the “FibroTrack Dist.zip”. Then, unzip the file by pressing extract all or extract here.





Following the extraction, open the folder and double click on the “FibroTrack Setup.exe”. Don’t delete “FibroTrack Setup-1a.bin” and “FibroTrack Setup-1b.bin”, these files are crucial for the software installation. Then, click on the Run anyway button.



Install the software by accepting the license.








The splash screen appears with a black console. It may take a few minutes to launch the software while importing all the modules for deep learning.

Welcome!

The app is initializing...
This may take a few minutes.



Please Wait...



Finally, a black console will appear with the following message.

```
Initializing the application...
please wait...


Welcome To

FibroTrack

Imports completed.
Please wait, you're almost there...
Launching the FibroTrack application...
```

Note: Don't close this console while you're using the software, otherwise you will close the software.

In the meanwhile, the application will be popped-up and cover the whole screen of your computer as follows:



Now, you are ready to use FibroTrack.

