


Guidelines for FibroTrack installation software for windows users

To install FibroTrack on windows, open


https://drive.google.com/file/d/1PvTxR_7k43wXrBx1fXlQomSxXSAPqRA2/view?usp=sharing

Click on the download button (Step1) and then Download anyway (Step 2). Finally, unzip the file by pressing “extract all” or “extract here” button (Step3).


Step 1



Step 2



Step 3



Step 4

Installation steps:


- Open the extracted “FibroTrack installation” folder and then open the “FibroTrack - Software Icons” folder.
- Find "FibroTrack Software.exe" (A)
- Double-click to run. You will have a pop up message: “Windows Protected your PC” click on “More info” (B).
- Click "Run anyway" (C).

Note: Do not delete any files from this folder - all contents are required for proper installation!


A

Name	Date modified	Type
▼ Today		
 FibroTrackSourceCodeFinal	06/08/2025 14:06	Application
 FibroTrackSourceCodeFinal.py	06/08/2025 14:06	PY File
 Group 1	06/08/2025 14:06	PNG File
 Group 2	06/08/2025 14:06	PNG File
 Odeh et al. Analyze Fibrosis	06/08/2025 14:06	JPG File
 Odeh et al. Define Pixel	06/08/2025 14:06	JPG File
 Odeh et al. Muscle Segmentation	06/08/2025 14:06	JPG File
 Splash Screen	06/08/2025 14:06	JPG File
 Splash_Screen_FibroTrack	06/08/2025 14:06	JPG File
 FibroTrack Software	06/08/2025 14:06	Application
 FibroTrack Software-1a.bin	06/08/2025 14:06	BIN File
 FibroTrackIcon	06/08/2025 14:06	JPG File
 best.pt	06/08/2025 14:06	PT File
 FibroTrack ICON	06/08/2025 14:06	ICO File

B






C




Step 5

Install the software by accepting the license.





The splash screen appears with a black console. It may take a few minutes to launch the software while importing all the modules for deep learning.



Finally, a black console will appear with the following message.


```
Initializing the application...
please wait...

Welcome To
FIBROTRACK

Imports completed.
Please wait, you're almost there...
Launching the FibroTrack application...
```

Note: Don't close this console while you're using the software, otherwise you will close the software. This process may take a few minutes.

In the meanwhile, the application will be popped up and cover the whole screen of your computer as follows:



Now, you are ready to use FibroTrack.